

- PCR1\_Forward (50 bp): 5'-TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCTACGGGNGGCWGCAG-3'
- PCR1\_Reverse (55 bp): 5'-GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGGACTACHVGGGTATCTAATCC-3'

What multiplexing capabilities are available with the protocol?

Nextera XT Index Kits use a dual-index strategy involving 12 (i7) and 8 (i5) indexes, for a total of 20 indexes. Together, they can generate 96 (12 × 8) different index combinations. For higher levels of multiplexing, Nextera XT v2 Index Kits allow up to 384 samples to be pooled in one flow cell lane. Nextera XT v2 Index Kits are available in four sets of 96 index combinations each, for a total of up to 384 samples per set.

How many samples can I pool per run?

The number of samples that can be pooled per run depends on the scope and goal of a given sequencing project. Generally, 15,000–100,000 reads per sample are sufficient for classification. For information about pooling strategies for a few samples, refer to the Nextera Low Plex Pooling Guidelines Technical Note<sup>6</sup>. Illumina recommends creating a sample sheet in Illumina Experiment Manager (IEM) before preparing libraries to confirm that the appropriate index combinations are used. If using MiSeq Control Software (MCS) 2.4.1 or 2.5 with real-time analysis (RTA) software 1.18.54, index confirmation is not necessary.

Sequencing

What are the recommended sequencing guidelines?

Illumina recommends using MiSeq<sup>®</sup> v3 sequencing chemistry. MiSeq v2 reagents are not recommended for 16S metagenomics studies, as they will affect quality scores and other metrics. A 5% PhiX spike-in is required with MiSeq v3 chemistry. Table 1 shows the expected yield for 16S sequencing runs using MiSeq v3 chemistry. With MiSeq Reporter software (2.3 or later), 50 bp of read overlap is required to take advantage of the read-stitching feature.

Table 1: MiSeq v3 Chemistry Specifications

Parameter	Value
Paired-end reads	44–50 M
Output	13.2–15 Gb
Quality scores	> 70% of bases higher than Q30 at 2 × 300 bp

Why are Q-scores lower in low-diversity libraries? How can I judge run performance?

Sequencing shotgun libraries generates intensity profiles, and the clusters on a flow cell represent roughly even proportions of A, C, G, and T at every cycle. In low-diversity libraries, such as those used for amplicon or 16S sequencing, the proportion of bases is unbalanced, so a large percentage of the clusters show the same base during each cycle. The high signals caused by the imbalance affect the base calling algorithm, resulting in low Q-scores even though the base calling accuracy is not necessarily poor. The 5% PhiX spike-in enables error rate calculations that allow verification of base calling accuracy over the course of the run, for all PhiX clusters, which can be extrapolated to the samples.

Does 16S sequencing require separate sequencing primers?

The required sequencing primers for the Illumina 16S metagenomics workflow are included in the MiSeq reagent cartridge and in HiSeq<sup>®</sup> 2500 Rapid-Run reagents. If using the HiSeq 2000 or HiSeq 2500 System for 16S sequencing, the TruSeq<sup>®</sup> Dual Index Sequencing Primer Box is required.

Is there a minimum percentage of sequence identity required for genus- and species-level classifications using the Illumina 16S metagenomics workflow?

The percentage of sequence identity needed depends on the organism under study. For more divergent species, longer reads can aid assembly and classification.



