

Part I. Welcome to Next-Generation Sequencing

NGS for *in vitro* Fertilization

Deciphering DNA sequences is essential for virtually all branches of biological research. With its unprecedented throughput, scalability, and speed, next-generation sequencing (NGS) enables more in-depth analysis of the genome than ever before, revealing new biological insights. In addition, NGS offers a faster, more comprehensive, streamlined method with more precise performance and higher throughput capabilities for obtaining DNA data than embryo morphology-based reviews. This new level of access to DNA information has catalyzed several important breakthroughs, making significant impacts in reproductive and genetic health studies, in particular in the field of *in vitro* fertilization (IVF).

As infertility rates rise and women wait until later in life to have children, so does the demand for IVF and assisted pregnancy. Unfortunately, IVF success rates remain low with only 30–35% of cycles resulting in a live birth.¹ Chromosome aneuploidy (abnormal number of chromosomes) is a major cause of IVF failure. Most embryos with aneuploidy do not implant, and those embryos that do often miscarry during the first trimester of pregnancy.^{2,3} Using preimplantation genetic screening (PGS) to identify euploid embryos, those embryos with a normal number of chromosomes, increases the chance that a viable embryo is selected for transfer during IVF, resulting in:

- Improved implantation rates^{4,5}
- Reduced spontaneous abortion^{4,6}
- Increased ongoing pregnancy rates^{4,6}

Using NGS to perform PGS delivers highly detailed, accurate, and scalable results for improved selection of euploid embryos over visual assessment alone.⁷

According to Dagan Wells, Associate Professor at the University of Oxford and Director at Reprogenetics-UK, “With next-generation sequencing technologies and PGS from Illumina, we can select embryos that are more likely to result in successful pregnancies”.

Simon Fishel, Managing Director of CARE Fertility, adds “PGS is a game changer for IVF clinics and labs. The ability to assess all 24 chromosomes in an embryo and the rapid delivery of the VeriSeq™ PGS Kit* empowers clinicians with accurate chromosomal insights into the health of the embryo to optimize IVF success.”

With NGS, researchers and clinicians have unprecedented access to genetic information. By enabling simultaneous analysis of multiple samples in fully automated workflows, NGS opens opportunities to increase efficiency, reduce costs, and shorten workflows. By generating highly specific and sensitive results, Illumina NGS solutions are setting new standards of success in the IVF landscape, enabling informed reproductive choices that can change patients’ lives.

Part II. What is Next-Generation Sequencing?

The Basic NGS Workflow


All Illumina NGS workflows include 4 basic steps.

1. Sequencing library preparation begins by creating short DNA or cDNA fragments with 5’ and 3’ adapters ligated (Figure 1A).
2. During cluster generation, the library is attached to an oligonucleotide lawn on the surface of a flow cell. Through bridge amplification, each library fragment acts as a seed to generate a clonal cluster containing thousands of identical fragments. Across the entire flow cell, millions to billions of clusters are formed (Figure 1B). Next, the templates are ready for sequencing by synthesis (SBS). SBS technology uses a proprietary reversible terminator–based method that detects single bases as they are incorporated into DNA template strands.⁸ Because all 4 reversible, terminator-bound dNTPs are present during each sequencing cycle, natural competition minimizes incorporation bias and greatly reduces raw error rates compared to other technologies.^{9,10}
3. The result is highly accurate base-by-base sequencing that virtually eliminates sequence-context-specific errors, even within repetitive regions and homopolymers (Figure 1C).
4. The newly identified sequence reads are then exported to an output file and aligned to a reference genome by sequencing alignment software (Figure 1D).

A detailed animation of SBS sequencing is available at www.youtube.com/watch?v=HMycqWhwB8E&list=UUxWMU29FF4kIG8YmQf6Zv0g


* The Illumina VeriSeq PGS Kit uses NGS to provide comprehensive, accurate screening of all 24 chromosomes simultaneously with an improved assay workflow and higher throughput compared to other methods.⁷

A. Library Preparation




NGS library is prepared by fragmenting a gDNA sample and ligating specialized adapters to both fragment ends.

B. Cluster Amplification



Library is loaded into a flow cell and the fragments hybridize to the flow cell surface. Each bound fragment is amplified into a clonal cluster through bridge amplification.


C. Sequencing



Sequencing reagents, including fluorescently labeled nucleotides, are added to the flow cell and the first base is incorporated. The flow cell is imaged and the emission wavelength and intensity are used to identify the base. This cycle is repeated "n" times to create a read length of "n" bases.

Figure 1: Next-Generation Sequencing (NGS) Workflow—The Illumina NGS workflow follows 4 basic steps: library preparation, cluster generation, sequencing, and data alignment.

D. Alignment & Data Analysis



Reads are aligned to a reference sequence with bioinformatics software. After alignment, differences between the reference genome and the newly sequenced reads can be identified.

Multiplexing

In addition to the rise of data output per run, the sample throughput per run in NGS has also increased over time. Multiplexing allows large numbers or batches of samples to be pooled and sequenced simultaneously during a single sequencing run (Figure 2). With multiplexed samples, unique index sequences are added to each DNA fragment during library preparation so that each read can be identified and sorted before final data analysis. Use of these indexes dramatically reduces the time-to-data for multi-sample studies and enables researchers to go from experiment to answer faster and easier than ever before.




Figure 2: Library Multiplexing Overview—A. Two distinct samples are attached to unique index sequences. Index sequences are attached during library preparation. B. Libraries for each sample are pooled together and loaded into the same flow cell lane. C. Samples are sequenced together during a single instrument run. All sequences are exported to a single output file. D. A demultiplexing algorithm sorts the reads into different files according to their indexes. E. Each set of reads is aligned to the appropriate reference sequence.

Part III. NGS Methods

NGS platforms enable a wide variety of applications, allowing researchers to investigate the genome, transcriptome, or epigenome of any organism. After sequencing libraries are prepared, the actual sequencing process remains fundamentally the same regardless of the application. Streamlined workflows, higher throughput capabilities, and customizable assays support easy portfolio expansion for any laboratory.

Sequencing for PGS

PGS performed using NGS follows a method known as targeted sequencing. With targeted sequencing, a preselected subset of genes is isolated and sequenced. These customized content sets are referred to as sequencing panels. Due to the tightly defined content, sequencing panels allow the genes within each sample to be sequenced to a greater depth than large-scale approaches like whole-genome or exome sequencing. This multigene approach confers several advantages:

- Decreases time-to-answer
- Minimizes issues with limited material for sequential testing
- More accurate and reliable results

Additional Sequencing Methods

The same platform used for sequencing in PGS can be used for various additional methods, including exome and whole-genome sequencing. This expands the utility of the system and enables users to expand their laboratory offerings.

Exome Sequencing

The exome, the protein-coding portion of the human genome, represents less than 2% of the genetic code, but contains ~85% of known disease-related variants.¹¹ With exome sequencing, the exonic regions of the genome are isolated and sequenced allowing researchers to focus on specific genetic regions of interest.

Whole-Genome Sequencing

Whole-genome sequencing (WGS) involves sequencing the full genome of a given organism. WGS can provide insight into the intronic (noncoding) regions of the genome that exome sequencing does not capture.

Part IV. Summary

Over the last decade, advances in NGS technology have led to an improved understanding of genomics, which, in turn, has led to new approaches to genetic screening. NGS has enabled significant strides in reproductive and genetic health, notably in IVF, setting a new standard in PGS, improving IVF success rates, and changing patients' lives.

Illumina is committed to providing the highest-quality data in the industry, setting a new standard of IVF success.

- Advancing trusted PGS and PGD genomic solutions
- Delivering accurate, timely, and reliable answers
- Offering educational resources and support

Learn more about Illumina NGS solution for reproductive and genetic health at www.illumina.com/NGSforIVF.

