# NovaSeq<sup>™</sup> 6000Dx Consumables



## Package Insert

FOR IN VITRO DIAGNOSTIC USE

## Intended Use

The Illumina NovaSeq 6000Dx Reagent v1.5 Kit (300 cycles)(S2 or S4) and the Illumina NovaSeq 6000Dx Buffer Cartridge (S2 or S4) are reagents and consumables intended for sequencing of sample libraries when used with validated assays. The kit is intended for use with the NovaSeq 6000Dx instrument and analytical software.

# **Principles of Procedure**

The Illumina NovaSeq 6000Dx S2 Reagent v1.5 Kit (300 cycles) and NovaSeq 6000Dx S4 Reagent v1.5 Kit (300 cycles) are single-use sets of reagents and consumables for sequencing on the NovaSeq 6000Dx Instrument. The Illumina NovaSeq 6000Dx S2 Reagent v1.5 Kit (300 cycles) and NovaSeq 6000Dx S4 Reagent v1.5 Kit (300 cycles) use libraries generated using validated assays. Sample libraries are captured on a flow cell and sequenced on the instrument using sequencing by synthesis (SBS) chemistry. SBS chemistry uses a reversible-terminator method to detect fluorescently labeled single nucleotide bases as they are incorporated into growing DNA strands. The number of sample libraries depends on the multiplexing supported by the upstream library preparation method.

The S2 Buffer Cartridge and the S4 Buffer Cartridge are single-use sets of buffers and wash solution that are used on the instrument during the sequencing run.

The NovaSeq 6000Dx Instrument Package Insert provides instructions for performing sequencing on the NovaSeq 6000Dx.

## Limitations of the Procedure

For in vitro diagnostic use only.

# **Product Components**

Table 1 NovaSeq 6000Dx Consumables

Kit Name	Illumina Catalog Number
NovaSeq 6000Dx S2 Reagent v1.5 Kit (300 cycles)	20046931
NovaSeq 6000Dx S4 Reagent v1.5 Kit (300 cycles)	20046933
NovaSeq 6000Dx S2 Buffer Cartridge	20062292
NovaSeq 6000Dx S4 Buffer Cartridge	20062293

### **Consumables**

Table 2 NovaSeq 6000Dx S2 Reagent v1.5 Kit (300 cycles)

Consumable	Quantity	Description
NovaSeq 6000Dx S2 Flow Cell	1	Single-use, paired-end glass flow cell
NovaSeq 6000Dx S2 Cluster Cartridge	1	Clustering reagents
NovaSeq 6000Dx S2 SBS Cartridge	1	Sequencing reagents

Table 3 NovaSeq 6000Dx S4 Reagent v1.5 Kit (300 cycles)

Consumable	Quantity	Description
NovaSeq 6000Dx S4 Flow Cell	1	Single-use, paired-end glass flow cell
NovaSeq 6000Dx S4 Cluster Cartridge	1	Clustering reagents
NovaSeq 6000Dx S4 SBS Cartridge	1	Sequencing reagents

Table 4 NovaSeq 6000Dx S2 Buffer Cartridge

Consumable	Quantity	Description
NovaSeq 6000Dx S2 Buffer cartridge	1	Buffers and wash solution

Table 5 NovaSeq 6000Dx S4 Buffer Cartridge

Consumable	Quantity	Description
NovaSeq 6000Dx S4 Buffer cartridge	1	Buffers and wash solution

# Storage and Handling

- 1. Storage conditions are defined as follows.
  - Room temperature is 15°C to 30°C
  - Refrigerated is 2°C to 8°C
  - Frozen is -25°C to -15°C
- 2. Reagent kit components are shipped under temperature control and are stable when stored at the indicated storage temperatures until the expiration date on the label.
- 3. The NovaSeq 6000Dx S2 Reagent v1.5 Kit (300 cycles) and NovaSeq 6000Dx S4 Reagent v1.5 Kit (300 cycles) reagent cartridges are stable for up to 4 hours when thawed in a controlled 19°C to 25°C water bath following specified procedure in the NovaSeq 6000Dx Instrument Package Insert.



4. Changes in the physical appearance of the reagents can indicate deterioration of the materials. If changes in the physical appearance occur, such as obvious changes in reagent color, or cloudiness apparent with microbial contamination, do not use the reagents.

## **Equipment and Materials Required, Sold Separately**

- NovaSeq 6000Dx, Illumina Catalog # 20068232
- NovaSeq 6000Dx Library Tube, Illumina Catalog # 20062290 (single) or Illumina Catalog # 20062291 (24 pack)

### **Lot Numbers**

Each reagent kit has a single lot number, which is referred to as the reagent kit lot number. Each box in the reagent kit is printed with the reagent kit lot number. Reagent kit components that are inside the boxes are printed with component-specific lot numbers that are different from the reagent kit lot number. Keep sequencing consumables stored in their boxes until ready for use to maintain kit lot association. See the Certificate of Analysis of the reagent kit for details about reagent part numbers and lot numbers.

# Warnings and Precautions



#### **CAUTION**

Federal law restricts this device to sale by or on the order of a physician or other practitioner licensed by the law of the State in which they practice, to use or order the use of the device.

- Some components of reagents provided by Illumina for use with the NovaSeq 6000Dx Instrument
  contain potentially hazardous chemicals. Personal injury can occur through inhalation, ingestion, skin
  contact, and eye contact. Wear protective equipment, including eye protection, gloves, and laboratory
  coat appropriate for risk of exposure. Handle used reagents as chemical waste and discard in
  accordance with applicable regional, national, and local laws and regulations. For additional
  environmental, health, and safety information, refer to the Safety Data Sheets (SDS) at
  support.illumina.com/sds.html.
- Failure to follow the procedures as outlined can result in erroneous results or significant reduction in sample quality.
- Use routine laboratory precautions. Do not pipette by mouth. Do not eat, drink, or smoke in designated work areas. Wear disposable gloves and laboratory coats when handling specimens and kit reagents. Wash hands thoroughly after handling specimens and kit reagents.
- Proper laboratory practices and good laboratory hygiene are required to prevent PCR products from contaminating reagents, instrumentation, and genomic DNA samples. PCR contamination can cause inaccurate and unreliable results.

## NovaSeq 6000Dx Consumables Package Insert

- To prevent contamination, make sure that pre-amplification and post-amplification areas have dedicated equipment and consumables (eg, pipettes, pipette tips, heat blocks, vortexers, and centrifuges).
- Index to sample pairing requires matching the index plate layout exactly. The DNA Prep with Enrichment
  Application automatically populates the index primers associated with the sample names, when entered
  during run setup. The user is advised to verify the index primers associate with samples before starting the
  sequencing run. Mismatches between the sample and plate layout results in loss of positive sample
  identification and incorrect result reporting.
- Do not operate the NovaSeq 6000Dx with any of the panels removed. Operating the instrument with any of the panels removed creates potential exposure to line voltage and DC voltages.
- Do not touch the flow cell stage in the flow cell compartment. The heater in this compartment operates between 22°C and 95°C and can result in burns.
- The instrument weighs approximately 1059 lbs. and could cause serious injury if dropped or mishandled.

## Instructions for Use

Refer to the NovaSeq 6000Dx Instrument Package Insert and applicable user documentation.

# **Performance Characteristics**

Refer to the NovaSeq 6000Dx Instrument Package Insert.

# **Revision History**

Document	Date	Description of Change
Document # 200025483 v00	August 2022	Initial release.



### **Patents and Trademarks**

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2022 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

### **Contact Information**



Illumina
5200 Illumina Way
San Diego, California 92122 U.S.A.
+1.800.809.ILMN (4566)
+1.858.202.4566 (outside North America)
techsupport@illumina.com





EC REP

Illumina Netherlands B.V. Steenoven 19 5626 DK Eindhoven The Netherlands

#### **Australian Sponsor**

Illumina Australia Pty Ltd Nursing Association Building Level 3, 535 Elizabeth Street Melbourne, VIC 3000 Australia

### **Product Labeling**

www.illumina.com

For a complete reference of symbols that appear on product packaging and labeling, refer to the symbol key at support.illumina.com on the *Documentation* tab for your kit.