

Streamline somatic variant interpretation

Identify insights faster while
scaling NGS operations for
somatic oncology research
applications





Move precision medicine forward

Streamlined

Accelerate time to research report by harnessing the power of automation, powerful visualization, and > 55 knowledge sources

Integrated

Remove unnecessary touchpoints and manual data movement by connecting data analysis with upstream sequencing to simplify and secure the comprehensive workflow

Powered for growth

Scale the volume and throughput of your NGS workflows without increasing headcount to keep pace with evolving science, technology, and demand

Connected Insights-Research harnesses > 55 knowledge sources to bring insights to diverse applications

Comprehensive knowledge in one place

- Connect LIMS, storage, pipelines, and third-party knowledge bases through APIs
- Integrate variant classifications, relevant drug labels, and pertinent clinical trials within a single view

Evidence-based integrated knowledge

- Leverage comprehensive solid and hematological cancer content, somatic gene variant annotations, and related content from multiple integrated knowledge bases:
 - OncoKB, a precision medicine knowledge base that includes FDA-recognized content from Memorial Sloan Kettering
 - JAX-CKB powered by The Jackson Laboratory, a trusted knowledge base with over 100,000 users worldwide

Flexible regionality

- Incorporate relevant lab- and region-specific practices within the user-defined workflow
- Support customized variant interpretation and reporting based on regional-tiering guidelines

Regular updates

- Access up-to-date, relevant variant annotations and insights, guidelines, and more with regular, automated content updates (as frequently as monthly for some knowledge sources)

Unlock meaning for a wide range of applications



Solid tumor testing (tissue)

Access comprehensive gene panels across multiple variant types, including TMB, MSI, and GIS to assess HRD



Liquid biopsy

Detect and analyze cancer variants with high analytical sensitivity and specificity using low levels of ctDNA



Hematological cancer

Streamline insight generation for myeloid leukemia, lymphoma, and other hematologic malignancies

Powerful features to streamline integration and adoption of a single-vendor workflow

Enterprise-level security and privacy standards

Protect the privacy of your genomic data with industry-leading global and local security standards

User-defined SOPs and advanced filters

Implement SOPs by creating sets of predefined custom filters or using the advanced filtering system

Streamlined workflow and collaboration

Facilitate teamwork, save time and remove manual data movement steps with auto-ingestion of variant data and autolaunch of Connected Insights-Research

Automated oncogenicity classification

Facilitate standardization within and across laboratories with automated, guideline-based oncogenicity classification, including a transparent evidence map

Broad portfolio of tests and variant types

Analyze comprehensive panels, exomes, genomes, or transcriptomes across a range of DNA and RNA variant and biomarker types, including TMB, MSI, and GIS to assess HRD

Lab-specific curation

Maintain a private knowledge base of your organization's curated data, including information about past variant interpretations and reporting

Automated custom research reporting

Customize, edit, and automatically populate draft research reports with minimal manual interaction

Powerful visualizations

Generate gene- and exon-level visualizations for DNA and RNA, variant QC, genome plots for structural variants, CNVs, B-allele ratio, fusion plots, and more

Integrate and streamline your workflows from library prep, sequencing, and data analysis

Learn more about the Connected Insights-Research at
illumina.com/connected-insights

White paper



GenomeWeb KOL white paper

Learn the thoughts of key opinion leaders on the current challenges and promise of NGS interpretation and reporting in clinical oncology



Data sheet



Connected Insights-Research data sheet

Read how Connected Insights-Research streamlines, integrates, and powers laboratories for scale and growth



Security brief



Connected Insights-Research security and compliance brief

Understand how extensive privacy, security, and compliance features make Connected Insights-Research a secure environment for sensitive data



illumina[®]

1.800.809.4566 toll-free (US) | +1.858.202.4566 tel
techsupport@illumina.com
www.illumina.com

© 2024 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.



Abbreviations

API	application programming interface	KOL	key opinion leader	SOP	standard operating procedure
CNV	copy number variant	LIMS	laboratory information management system	SV	structural variant
ctDNA	circulating tumor DNA	MSI	microsatellite instability	TMB	tumor mutational burden
GIS	genomic instability score	NGS	next-generation sequencing	WES	whole-exome sequencing
HRD	homologous recombination deficiency	SNV	single nucleotide variant	WGS	whole-genome sequencing
Indels	insertion/deletion				